


Device alternatives for interventions in pediatric cardiology and congenital heart disease

The following file contains information about different devices and especially device alternatives for interventional procedures in pediatric cardiology and congenital heart disease. The information is continuously collected from members of the Interventional WG of the AEPC following an initiative of the Council of the WG. It represents personal experiences of the users and does not contain any recommendations of the interventional WG of the AEPC. The interventional WG of the AEPC does not take any responsibility about the information in this file.

- Stents
- Sheaths
- Catheters
- Balloons
- Devices
- Others



Stents

Bare metal coronary stents

- [Avantgarde](#)


Brands

- Alternatives to Formula:
 - [Palmaz blue](#)
 - [Palmaz Genesis](#)
 - [VBX Balloon expandable](#)
 - [Besmooth](#)
 - [RX Herculink Elite](#)
- Alternatives to ev3:
 - [Optimus](#)
 - [Numed CP](#)
 - [Bentley Graft aortic](#)

Sizes

- [2,25-4,5](#)
- [4-10 mm](#)
- [6-18 mm](#)
- [10-24 mm](#)


Bare metal coronary stents not drug eluting



- Avantgarde
 - [Alvimedica](#)
 - 2,25-4,5 mm diameter
 - 8-31 length
 - 5F GC, 0,014 guidewire

t.krasemann@erasmusmc.nl

Premounted stents – 4-10mm alternatives to Cook Formula Stents



- Palmaz Blue
 - [Cordis](#)
 - 4-7mm diameter
 - 10-24 length
 - 5F sheath, 0,014 or 0,018 guidewire
 - Can be overdilated more than mentioned in the brochure
 - fractures at 10mm

t.krasemann@erasmusmc.nl
georgiev@dhm.mhn.de

- Palmaz Genesis
 - [Cordis](#)
 - diameter: 5-10 mm
 - Length: 17-80 mm
 - Sheath: 6-7 F
 - Recommended inflation pressure : 8 atm


gianfranco.butera@opbg.net

- VBX Balloon Expandable
 - [Gore](#)
 - Covered stents
 - Size: 5-11 mm
 - Length: 15-79 mm
 - Post-dilated diameter ~ +5 mm
 - Sheath max 8 Fr

gianfranco.butera@opbg.net

Premounted stents – 4-10mm

alternatives to Cook Formula Stents




- Besmooth
- [Bentley InnoMed Hechigen](#)
- Diameter 5-10 mm
- Length 18-56 mm
- Premounted
- Flexible
- Low profile 6 Fr
- Failure to post dilate beyond 12mm

guiti.milani@gmail.com

- RX Herculink Elite
 - [Abbott](#)
- Diameter 4-7 mm
- Length 12-15-18mm
- Co Cr stent
- 6 Fr sheath
- ? redilatable

guiti.milani@gmail.com


Balloon expandable stents – 6-30mm alternatives to EV3 MegaLD and MaxLD



Optimus L, XL, XXL

- [Andratec](#)
- Dilatable from 6 up to 20mm (L), 24 (XL) and 30 (XXL)mm
- Sheath=balloon+1F
- Different lengths
- Covered version available as well

georgiev@dhm.mhn.de



Balloon expandable stents – 6-30mm

Numed


- CP stent 8 Zig
 - Size: 12-24 mm,
 - Stent length: 16-60 mm
- CP stent 10 zig
 - Size: 26-30 mm
 - Stent length: 39-60 mm
- Sheath = balloon +2F
- Considerable capacity for expansion
- Covered version available

gianfranco.butera@opbg.net

BeGraft aortic

- Bentley
- Premounted covered stent (ePTFE)
- Size: 12-24mm
 - At nominal pressure 12, 14 mm->max. diameter 20 mm;
 - At nominal pressure 16, 18 mm-> max. diameter 24 mm
 - At nominal pressure 20, 22, 24 mm-> max. diameter 30 mm
- Length: 19-59 mm
- Sheath: Ø 12 mm -> 9F; Ø 14, 16 mm -> 11 F; Ø 18, 20, 22, 24 mm-> 14 F


gianfranco.butera@opbg.net



Long sheaths 6-26F

Models

- [Dryseal](#)
- [Flexor](#)
- [Destination](#)
- [D`Vill](#)
- [APT](#)
- [Fortress](#)
- [Occlutech](#)



Long sheaths 6-26F

Dryseal

- [Gore](#)
- Sizes 10-26F, increment 2F
- Length 33 and 65 cm
- Good manoeuvrability
- Minimal shoulder between dilator and sheath
- Requires 0.035 guidewire

t.krasemann@erasmusmc.nl

• Flexor Check-Flo[®] Valve


- [Cook](#)
- Size: 4-12 Fr
- Length: 45-55-70-90-110 mm.
- Accepts wire: 0.018/0.035
- Tip Configuration: straight/multipurpose
- Good manoeuvrability

gianfranco.butera@opbg.net

• Destination[™] Guiding Sheath

- [Terumo](#)
- Fr: 5-8 Fr
- Length: 45-65-90 mm.
- Maximum guidewire size: 0.035
- Reinforcement inside the sheath with stainless steel coil
- Tip Configuration:
 - Straight
 - Hockey Stick
 - Multipurpose
 - RDC
 - LIMA

gianfranco.butera@opbg.net



Long sheaths 6-26F

D'VILL™ Introducer

- [Brown](#)
- Long length
- Size: 12-14Fr
- Length: 30cm, 65cm, and 85cm
- Good manoeuvrability

gianfranco.butera@opbg.net

APT long sheath


- [APT Medical](#)
- Size: 4-20Fr

guiti.milani@gmail.com

Fortress

- [Biotronic](#)
- Size: 4-8Fr
- Length: 45, 65, 90 cm
- Reinforced
- Relatively stiff dilator
- Very good manoeuvrability, especially over a balloon

Joerg.Michel@med.uni-tuebingen.de




Long sheaths 6-26F

Occlutech sheaths

- [Occlutech](#)
- Sheaths for device implantation, very useful for stent implantation (off label)
- Size: 6-14Fr
- Good manoeuvrability


georgiev@dhm.mhn.de

Catheters



Types

- [Wedge](#)



Wedge catheters


alternatives to Arrow Wedge Catheter

Corodyn P1

- [Brown](#)
- Wedge Pressure
- Balloon Flotation Catheter
- Single Pressure Lumen
- Sizes 5, 6, and 7Fr

gianfrancobutera@libero.it

Balloons




Models

- [Numed Rashkind balloons](#)
- [Edwards Rashkind balloons](#)
- [Tashak](#)
- [VACS II](#)
- [Z-Med](#)
- [BiB](#)
- [VACS III](#)
- [AltoSa](#)
- [Atlas Gold](#)
- [Dorado](#)
- [Conquest](#)
- [Athletis](#)
- [True Dilatation](#)
- [Cristal](#)
- [Advance LP](#)
- [Mustang](#)
- [Wolverine](#)
- [Angiosculpt](#)
- [Agent](#)
- [In.pact Admiral](#)
- [Armada](#)
- [Euca PW](#)

Types

- [Rashkind Balloons](#)
- [Valve dilatations](#)
- [Stent implantation](#)
- [High pressure](#)
- [Ultrahigh pressure](#)
 - <12mm
 - >12mm
- [Cutting balloons](#)
- [Drug coated balloons](#)

Rashkind Balloons




- NuMed Z5 and Z6
- Over the wire, but can be used without
- 9,5mm and 13,5mm
- Non compliant
- 6F sheath
- Balloon deflation sometimes incomplete-cave sheath!

t.krasemann@erasmusmc.nl

- Edwards
- No wire lumen
- Non compliant
- 6F sheath
- Gut steerable and visible under echo
- Exceptional approval in Germany

georgiev@dhm.mhn.de

Balloons for valve dilatations




- Tyshak
 - [Numed](#)
 - 2-30mm diameter
 - 1-10cm length
 - Low profiles
 - Excellent trackability
 - Low pressures

georgiev@dhm.mhn.de

- VACS II
 - [Osypka](#)
 - 4-30mm diameter
 - 2-6cm length
 - Low profiles
 - Very good trackability
 - Low pressures

georgiev@dhm.mhn.de

Balloons for stent implantation



- Z-med
 - [Numed](#)
 - 2-40mm diameter
 - 2-8cm length
 - Profile starting from 6F

georgiev@dhm.mhn.de


- BiB
 - [Numed](#)
 - Balloon-in-balloon
 - 12-30mm diameter
 - 2,5-6,5cm length
 - Profile starting from 8F
 - Inner balloon 50% from outer

georgiev@dhm.mhn.de

- VACS III
 - [Osypka](#)
 - 5-30mm diameter
 - 2-6cm length
 - High pressure dilatations possible

georgiev@dhm.mhn.de

Balloons for stent implantation




- AltoSa
 - [Andratec](#)
 - 12-30mm diameter
 - 2-6cm length
 - Low profile starting from 7F

georgiev@dhm.mhn.de

- Gemini AltoSa
 - [Andratec](#)
 - Balloon-in-balloon
 - 14-30mm diameter
 - 3-6cm length
 - Profile 8-12F
 - Inner balloon about 60% from outer

georgiev@dhm.mhn.de

High- and ultrahigh pressure balloons



Atlas Gold

- [Becton Dickinson](#)
- Kevlar balloon
- 12-26mm
- Up to 18 Bar (off label more, ruptures punctually)

georgiev@dhm.mhn.de

Dorado

- [Becton Dickinson](#)
- Wooven balloon
- 3-10mm


georgiev@dhm.mhn.de

Conquest

- [Becton Dickinson](#)
- Wooven balloon
- 6-10mm
- Up to 40 bar

georgiev@dhm.mhn.de

High- and ultrahigh pressure balloons



Athletis

- [Boston Scientific](#)
- Kevlar balloon
- 4-12mm
- Up to 40 Bar

georgiev@dhm.mhn.de

True Dilatation

- [Becton Dickinson](#)
- 18-28mm

georgiev@dhm.mhn.de

Cristal Balloon

- [Balt](#)
- 12-28mm
- Up to 10 bar
- Semi-compliant
- Low profile
- Size: 8-20 mm
- Length: 20-40 mm
- Sheath: 6-9 Fr

gianfranco.butera@opbg.net

High pressure balloons 4-12mm



Advance 14LP

- 0.014, 2-4mm/2-20cm, 4 F

Advance 18LP

- 0.018, 4-9mm/2-20 cm, 5-6 Fr

Advance 35LP

- 0.035, 3-12mm/2-20 cm, 5-7 Fr

- [Cook](#)

- Up to 14 bar

gianfranco.butera@opbg.net

Mustang balloon

- [Boston Scientific](#)
- Wire: 0.035
- Size: 3-12 mm
- Length: 2-10 cm
- Sheath: 5-7 F
- High pressure balloon (up to 24 atm)
- Good for stent implantation


gianfranco.butera@opbg.net

Armada balloon

- [Abbott](#)
- Diameter 3-14mm
- Length 20-250 mm
- 5-7 Fr sheath
- Catheter length 80-135 cm
- Nominal pressure 6-9 ATM
- Burst pressure 10-20 ATM

guiti.milani@gmail.com


High pressure balloons 4-12mm



Euca PW

- [Elaine Medical](#)
- Diameter 4-12 mm
- Length 20-200mm
- 5-7 Fr
- PTA balloon
- OTW 0.035

guiti.milani@gmail.com



Cutting balloons

Wolverine


- [Boston Scientific](#)
- Blade Length: 6/10/15 mm
- Size: 2-4 mm

gianfranco.butera@opbg.net

Angiosculpt

- [Philips](#)
- Size: 2-8 mm
- Length: 10-200 mm
- Sheath: 5/6F
- Wire: 0.014/0.018

gianfranco.butera@opbg.net




Drug coated balloons


- AGENT™, Boston Scientific
 - Size: 2-4 mm
 - Length: 12/15/20/30 mm
 - Wire: 0.014
- IN.PACT Admiral, Medtronic
 - Size: 4-7 mm
 - Length: 40,60,80,120,150,200, and 250 mm
 - Over-the-wire
 - Wire: 0,035

gianfranco.butera@opbg.net

Devices



- [ASD](#)
- [PDA](#)
- [VSD](#)



Device: ASD

Amplatzer ASD

- Waist: 4-40 mm
- Sheath: 6-12 Fr

Occlutech ASD

- Waist: 4-40 mm
- Sheath: 7-12 Fr


Cocoon ASD

- Waist: 8-42 mm
- Sheath: 6-14 Fr

Gore septal occluder (GCO)

- Size:
 - 27 mm (defect 8-15mm)
 - 32 mm (defect 13-20 mm)
 - 37 mm (defect 18-25 mm)
 - 44 mm (defect 23-30 mm)
 - 48 mm (defect 28-35 mm)

Sheath: 10-14 Fr



Device: PDA

Amplatzer duct occluder (ADO)

- Different size
- Standard length
- Sheath: 5-7 Fr

Occlutech PDA

- Different size
- Sheath: 6-9 Fr
- Shank: standard/long

PFM coil


- 4x4mm and 4x5 mm can be implanted from the aorta.

Piccolo occluder

- 3-5mm size
- 2-6mm length
- Sheath: 4 Fr
- Only approved devices for premies > 700g

Medtronic MVP

- 3-7mm size
- 12mm length
- Microcatheter
- Suitable for premature PDAs



Device: VSD


Amplatz-VSD-MUSC

- Waist: 4-18 mm
- Sheath: 6-9F

Konar MF VSD

- Size: 2-10 mm
- Sheath: 4-7Fr

Others



- Bioptomes
- Bailys system

Contacts to contributors

- Stanimir Georgiev, German Heart Center Munich georgiev@dhm.mhn.de
- Thomas Krasemann, Erasmus MC, Rotterdam t.krasemann@erasmusmc.nl
- Gianfranco Butera, Bambino Jesu, Rome gianfranco.butera@opbg.net
- Guiti Milani, Necker, Paris guiti.milani@gmail.com
- Jörg Michel, Tübingen, Germany Joerg.Michel@med.uni-tuebingen.de